



UNIVERSITÀ
DI TRENTO


Measurement of the curvature of a mirror using a Michelson interferometer

Candidate: Anna Do'

Supervisor: Giacomo Lamporesi


Bachelor's Degree in Physics

20 November 2024



The topics covered in this presentation are:

- Michelson interferometer
- Two theoretical models:
geometric optics and
interference of Gaussian
beam
- Data analysis
- Phase shift method




The electric field can be expressed in the complex notation as:


$$\vec{E} = \vec{E}_0 e^{i\Phi}$$

If \vec{E}_1 and \vec{E}_2 are the electric in the two arms, the total field is $\vec{E}_{tot} = \vec{E}_1 + \vec{E}_2$, then the intensity on the screen is proportional to:


$$I \propto \langle |\vec{E}_{tot}|^2 \rangle = 2I_0^2 [1 + \cos(\Phi_1 - \Phi_2)] = 2I_0^2 [1 + \cos(k\Delta z)]$$




Experimental set-up



Experimental set-up





The optical path difference in the two arms of the interferometer is:

$$\Delta z = 2d - (d + \overline{NM} + \overline{ML}) = \frac{2(1 - \cos \theta)[d + d \cos \theta + R \cos^2 \theta]}{2 \cos^2 \theta - 1} \simeq$$

$$\simeq \frac{r^2}{R^2}(2d + R) \simeq \frac{r^2}{R} \rightarrow I \propto \left[1 + \cos \left(\frac{kr^2}{R} \right) \right]$$

The equation of a Gaussian beam is:

$$E = E_0 \frac{w_0}{w(z)} \exp\left(\frac{-r^2}{w^2(z)}\right) \exp\left[ik\left(z + \frac{r^2}{2R(z)}\right) - i \tan^{-1}\left(\frac{z}{z_0}\right)\right]$$

where:

$$z_0 = \frac{\pi w_0^2}{\lambda} \quad w(z) = w_0 \sqrt{1 + \left(\frac{z}{z_0}\right)^2} \quad R(z) = z \left[1 + \left(\frac{z_0}{z}\right)^2\right]$$

Let us call Beam₁ the beam reflected by the spherical mirror and Beam₂ the beam coming from the planar mirror:

	w_0	z_0	$w(z)$	$R(z)$
Beam ₁	$\simeq 0.025 \text{ mm}$	$\simeq 3.8 \text{ mm}$	$w_{01} z_1 / z_{01}$	z_1
Beam ₂	$\simeq 5 \text{ mm}$	$\simeq 150 \text{ m}$	w_{02}	z_{02}^2 / z_2

The intensity of the electric magnetic field on the screen will be:

$$I \propto |E_{tot}|^2 =$$

$$= E_0^2 \left\{ \exp\left(-\frac{2r^2}{w_{02}^2}\right) + \left(\frac{z_{01}}{z_1}\right)^2 \exp\left(-\frac{2r^2 z_{01}^2}{w_{01}^2 z_1^2}\right) + \right. \\ \left. + \left(\frac{z_{01}}{z_1}\right) \exp\left[-r^2 \left(\frac{1}{w_{02}^2} + \left(\frac{z_{01}}{w_{01} z_1}\right)^2\right)\right] 2 \cos\left[k(z_2 - z_1) + \frac{kr^2}{2} \left(\frac{z_2}{z_{02}^2} + \frac{1}{z_1}\right)\right] \right\} \simeq$$

$$\simeq E_0^2 \left[\exp\left(-\frac{2r^2}{w_{02}^2}\right) + \left(\frac{z_{01}}{z_1}\right)^2 + 2 \left(\frac{z_{01}}{z_1}\right)^2 \exp\left(-\frac{r^2}{w_{02}^2}\right) \cos\left(\frac{kr^2}{R}\right) \right]$$


The fit function is:


$$I = A [B^2 + \exp(-2ar^2) + 2B \exp(-ar^2) \cos(br^2 + c)]$$

where from the previous equation
we have that:

$$a = \frac{1}{w_{02}^2} \quad b = \frac{2\pi}{\lambda(R + 2d)}$$

$$B = \frac{z_{01}}{z_1}$$





$$R = \frac{2\pi}{b\lambda} = (2.57 \pm 0.01) \text{ m}$$

Let us consider the intensity on the camera at different position of the planar mirror ($d + \delta\Phi_z$):

$$I = A(x, y) \left[1 + B(x, y) \cos \left(\Phi(x, y) + \delta\Phi_z \right) \right] + off(x, y)$$

The surface of the mirror $h(x, y)$ is related to the phase $\Phi(x, y)$:

$$h(x, y) = \frac{\lambda}{4\pi} \Phi(x, y)$$

For fixed value of r , we observe that the intensity is a sinusoidal function of $\delta\Phi_z$:

$$I(\delta\Phi_z) = A \sin(\delta\Phi_z + \Phi) + off = A \sin(nv\Delta t + \Phi) + off$$





Figure: The first image is for the pixel with coordinates (0,0) and the second for the pixel with coordinates (20,50).

Analysis of the data performed by the code:

- get Φ from the data, there are two option:
 - sinusoidal fit
 - fft : fast Fourier transform
- unwrap the phase Φ
- fit to get the parameter that describe the surface, there are two option:
 - $\left[A \left((y - y_c) \cos(\alpha) + (x - x_c) \sin(\alpha) \right)^2 + B \left((x - x_c) \cos(\alpha) - (y - y_c) \sin(\alpha) \right)^2 \right] + o$
 - Zernike polynomials